COLUMBIA DOCTORS TARRYTOWN

NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591 CD Submission 06/18/2021

OWNER	A
COLUMBIA UNIVERSITY MEDICAL	ARF

CENTER (CUMC) 155 WHITE PLAINS ROAD

TARRYTOWN, NY 10591 PHONE: 914.333.9800

ARCHITECT

RRAY ARCHITECTS 470 PARK AVENUE SOUTH, 11TH

New York, NY, 10016 PHONE: 212.689.3110

STRUCTURAL ENGINEER

REUTHER + BOWEN

DUNMORE, PA 18512 PHONE: 570.496.7020


326 WARD STREET

MEP ENGINEER

LORING CONSULTING ENGINEERS,

360 WEST 31ST STREET

NEW YORK, NY 10001 PHONE: 646.674.6100







NEW YORK, NY. 10016 212-689-3110

CONSULTANTS:

STRUCTURAL ENGINEER

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER


LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001 PHONE: 646.674.6100


OWNER:

COLUMBIA DOCTOR'S TARRYTOWN

PROJECT: NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591







SHEET TITLE: COVER SHEET & DRAWING LIST


SEAL:

DATE:06/18/2021 CON/REF No. CONTRACT No. SCALE: AS NOTED PROJECT No. 6109 CHECKED: CH DRAWN: KU


SHEET NO. T-000.00

DWG OF 47

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be prosecuted under the law.



TITLE T-000	.00	COVER SHEET & DRAWING LIST
GENERAL	.00	COVER SHEET & DIVAVVIIVO LIST
G-001	.00	CODE REQUIREMENTS, SYMBOLS LEGEND, FLOOD MAP DATA AND GENERAL NOTES
G-002	.00	UL DETAILS
G-003	.00	UL DETAILS
G-004	.00	UL DETAILS
G-010 G-011	.00	ICRA PLAN AND NOTES PROJECT PHASING AND MRI RIGGING
LIFE SAFET		TROJECT FRANKS AND WIRTRIGGING
LSC-101	.00	LIFE SAFETY & OCCUPANCY PLANS AND NOTES
DEMOLITIO		
AD-100 STRUCTUR	.00	DEMOLITION PLAN
S-001	.00	GENERAL STRUCTURAL NOTES AND DESIGN CRITERIA
S-100	.00	PARTIAL FOUNDATIONS PLAN AND DETAILS
S-101	.00	2ND FLOOR AND ROOF FRAMING PLAN AND DETAILS
ARCHITECT		ADCLUTECTUDAL DI ANI
A-100 A-300	.00	ARCHITECTURAL PLAN ROOF PLAN
A-301	.00	ROOF DETAILS
A-400	.00	EQUIPMENT PLAN
A-500	.00	REFLECTED CEILING PLAN & CEILING DETAILS
A-600 A-720	.00	INTERIOR ELEVATIONS PARTITION TYPES & RATED ASSEMBLY
A-720 A-721	.00	FRAMING DETAILS
A-722	.00	INTERIOR DETAILS
A-800	.00	FINISH PLAN AND SCHEDULE
A-900	.00	DOOR TYPES, FRAME DETAILS, SCHEDULE & BORROWED LIGHT DETAILS, SCHEDULE
VENDOR		DE IAILO, OOI ILDULL
V-101	.00	GE MRI-C1
V-102	.00	GE MR-C2
V-103	.00	GE MR-A1
V-104 V-105	.00	GE MR-A2 GE MR-A3
V-105 V-106	.00	GE MRI-A4
V-107	.00	GE MR-A5
V-108	.00	GE MRI-A6
V-109 V-110	.00	GE MRI-A7 GE MRI-A8
V-110 V-111	.00	GE MR-A9
V-112	.00	GE MRI-A10
V-113	.00	GE MRI-A11
V-114	.00	GE MR-S1
V-115 V-116	.00	GE MRI-S2 GE MRI-S3
V-117	.00	GE MRI-M1
V-118	.00	GE MR-M2
V-119	.00	GE MRI-M3
V-120 V-121	.00	GE MRI-M4 GE MRI-M5
V-122	.00	GE MR-E1
V-123	.00	GE MR-E2
V-124	.00	GE MRI-E3
V-125 V-126	.00	GE MR-E4 GE MRI-E5
V-127	.00	GE MR-E6
MECHANICA	٩L	
M-101	.00	DEMOLITION PART PLAN, HVAC PIPING AND DUCTWORK
M-201 M-202	.00	FIRST FLOOR HVAC PART PLAN, DUCTWORK 2ND FL AND ROOF HVAC PART PLAN PIPING AND CRYOGENIC VENT
M-001	.00	MECHANICAL SYMBOL LIST, NOTES, AND ABBREVIATIONS
ELECTRICA		
E-001	.00	ELECTRICAL SPECIFICATION SHEET
E-002 E-003	.00	ELECTRICAL SPECIFICATION SHEET ELECTRICAL NOTES
E-101	.00	LEVEL1 POWER DEMOLITION PART PLAN
E-102	.00	LEVEL1 LIGHTING DEMOLITION PART PLAN
E-201	.00	LEVEL1 POWER PART PLAN
E-202 E-301	.00	ROOF LEVEL POWER PART PLAN LEVEL 1 LIGHITNG PART PLAN
E-301 E-401	.00	ELECTRICAL RISER DIAGRAM
E-501	.00	ELECTRICAL SCHEDULES
E-601	.00	ELECTRICAL DETAILS
E-602 PLUMBING	.00	ELECTRICAL DETAILS
P-001	.00	PLUMBING NOTES, SYMBOLS, SCHEDULES & DETAILS
P-002	.00	PLUMBING SPECIFICATIONS
P-101	.00	PLUMBING DEMOLITION PLAN
P-201 FIRE ALARN	.00 <i>I</i>	PLUMBING PLANS
FA-001	.00	FIRE ALARM COVER SHEET
FA-101	.00	LEVEL1 FIRE ALARM DEMOLITION PART PLAN
FA-301	.00	LEVEL1 FIRE ALARM PART PLAN
FIRE PROTE		SDDINKI ED NOTES SYMBOLS SCHEDULES & DETAILS
SP-001 SP-002	.00	SPRINKLER NOTES, SYMBOLS, SCHEDULES & DETAILS SPRINKLER SPECIFICATIONS
		SPRINKLER DEMOLITION PLAN
SP-101	.00	OF MINKELY DEMOCITION FEAN



CONTRACTOR GENERAL NOTES



PHONE: 570.496.7020

MEP ENGINEER

DUNMORE, PA 18512


LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001 PHONE: 646.674.6100


OWNER:

COLUMBIA DOCTOR'S TARRYTOWN

PROJECT: **NEW MRI**

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591





DESCRIPTION REVISIONS/ISSUES

SHEET TITLE: CODE REQUIREMENTS, SYMBOLS LEGEND, FLOOD MAP DATA AND GENERAL NOTES

SEAL:


DATE:06/18/2021 CON/REF No. CONTRACT No. **SCALE:** AS NOTED **PROJECT No.** 6109 CHECKED: CH DRAWN: KU

SHEET NO.

G-001.00

DWG OF 47

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be



Architect of Record:

ARRAY-ARCHITECTS.COM

470 PARK AVE SOUTH, 11th FLOOR NEW YORK, NY. 10016 212-689-3110

STRUCTURAL ENGINEER

CONSULTANTS:

REUTHER + BOWEN
326 WARD STREET
DUNMORE, PA 18512
PHONE: 570.496.7020

MEP ENGINEER

LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001 PHONE: 646.674.6100


OWNER:

COLUMBIA DOCTOR'S TARRYTOWN

PROJECT:
NEW MRI

55 WHITE PLAINS ROAD

W



NO. DESCRIPTION DATE
REVISIONS/ISSUES

SHEET TITLE:
UL DETAILS

SEAL:


DATE:06/18/2021
CON/REF No.
CONTRACT No.
SCALE: AS NOTED
PROJECT No. 6109
CHECKED:CH
DRAWN: KU


SHEET NO.

G-002.00

This document is a copyright protected instrument of service,


property of Array Architects and licensed for use in the title







CD Submission 06/18/2021


This document is a copyright protected instrument of service property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be prosecuted under the law.



written permission of Array Architects is illegal and will be prosecuted under the law.








470 PARK AVE SOUTH, 11th FLOOR

W - WEST BUILDING E - EAST BUILDING


CON/REF No. CONTRACT No. SCALE: AS NOTED PROJECT No. 6109 CHECKED: CH DRAWN: KU


AD-100.00

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be prosecuted under the law.



written permission of Array Architects is illegal and will be prosecuted under the law.





Architect of Record:

GENERAL CONSTRUCTION NOTES

ALL CONSTRUCTION AND MATERIALS FOR MRI ROOM 2/ ZONE 4 - W189C AND MRI

ARRAY-ARCHITECTS.COM 470 PARK AVE SOUTH, 11th FLOOR

CONSULTANTS:

STRUCTURAL ENGINEER

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER


LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001 PHONE: 646.674.6100

OWNER:

COLUMBIA DOCTOR'S TARRYTOWN

NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591



KEY PLAN

DESCRIPTION DATE **REVISIONS/ISSUES**

SHEET TITLE: **ROOF DETAILS**

SEAL:

CON/REF No. CONTRACT No. SCALE: AS NOTED PROJECT No. 6109 CHECKED: CH DRAWN: KU

DATE:06/18/2021


A-301.00


This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without

written permission of Array Architects is illegal and will be

prosecuted under the law. C COPYRIGHT - 2020 ARRAY-ARCHITECTS P.C.

CD Submission 06/18/2021





Architect of Record:

ARRAY-ARCHITECTS.COM

CONSULTANTS:

STRUCTURAL ENGINEER

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020


MEP ENGINEER

LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001

COLUMBIA DOCTOR'S TARRYTOWN

NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591



DESCRIPTION REVISIONS/ISSUES

SHEET TITLE: REFLECTED CEILING PLAN & CEILING DETAILS

CON/REF No. CONTRACT No. **SCALE:** AS NOTED PROJECT No. 6109 CHECKED: CH DRAWN: KU


DATE:06/18/2021

A-500.00

C COPYRIGHT - 2020 ARRAY-ARCHITECTS P.C.

OF **47**

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be prosecuted under the law.



GENERAL CONSTRUCTION NOTES

ALL CONSTRUCTION AND MATERIALS FOR MRI ROOM 2/ ZONE 4 - W189C AND MRI CONTROL ROOM / ZONE 3 - W204F ARE TO BE NON-FERROMAGNETIC. NON-FERROMAGNETIC MATERIALS INCLUDING BUT NOT LIMITED TO REMOVABLE COVERS, SWITCH BOX COVER PLATES, LIGHT FIXTURE COMPONENTS, MOUNTING SCREWS, ETC.

TO MINIMIZE POTENTIAL HAZARD IN THE MAGNET ROOM.

MRI ZONE SAFETY NOTES

SUITES FOR MRI EQUIPMENT SHALL BE PLANNED TO CONFORM TO THE FOUR-ZONE SCREENING AND ACCESS CONTROL PROTOCOLS IDENTIFIED IN THE CURRENT EDITION OF THE AMERICAN COLLEGE OF RADIOLOGY'S "GUIDANCE DOCUMENT FOR SAFE MRI PRACTICES." AMERICAN COLLEGE OF RADIOLOGY (ACR) ZONING DESCRIPTIONS:

ZONE 1: ALL AREAS FREELY ACCESSIBLE TO THE GENERAL PUBLIC WITHOUT SUPERVISION. MAGNETIC FRINGE FIELDS IN THIS AREA ARE LESS THAN 5 GAUSS

ZONE 2: STILL A PUBLIC AREA, BUT THE INTERFACE BETWEEN UNREGULATED ZONE 1 AND THE STRICTLY CONTROLLED ZONES 3 AND 4. MR SAFETY SCREENING TYPICALLY OCCURS HERE UNDER TECHNOLOGIST SUPERVISION.

ZONE 3: AN AREA NEAR THE MAGNET ROOM WHERE THE FRINGE, GRADIENT, OR RF MAGNETIC FIELDS ARE SUFFICIENTLY STRONG TO PRESENT A PHYSICAL

HAZARD TO UNSCREENED PATIENTS AND PERSONELL. ZONE 4: SYNONYMOUS WITH THE MR MAGNET ROOM ITSELF. HAS THE HIGHEST FIELD (AND GREATEST RISK) AND FROM WHICH ALL FERROMAGNETIC OBJECTS

GENERAL NOTES

6. SEE SHETT A-011 FOR MRI RIGGING PATH AND SAFETY ZONE PLAN.

. SEE SHEET T-000 FOR COVER SHEET & DRAWING LIST.

2. SEE SHEET G-001 FOR CODE REQUIREMENTS, SYMBOLS LEGEND, FLOOD MAP DATA AND GENERAL NOTES.

3. SEE SHEET LSC-101 FOR LIFE SAFETY & NOTES.

4. SEE SHEET DM-100 FOR DEMOLITION PLAN. 5. SEE SHEET A-010 FOR ICRA PLAN AND NOTES.

MUST BE INCLUDED.

7. SEE SHEET A-100 FOR ARCHITECTURAL PLANS. 8. SEE SHEET A-300 FOR ROOF PLAN.

9. SEE SHEET A-400 FOR EQUIPMENT PLAN. 10. SEE SHEET A-500 FOR REFLECTED CEILING PLAN & CEILING DETAILS.

11. SEE SHEET A-721 FOR PARTITION TYPES & RATED ASSEMBLY.

12. SEE SHEET A-722 FOR INTERIOR DETAILS.

13. SEE SHEET A-800 FOR FINISH PLAN AND SCHEDULE.

14. SEE SHEET A-900 FOR DOOR TYPES, FRAME DETAILS, SCHEDULE & BORROWED LIGHT DETAILS, SCHEDULE.

ARCHITECTURAL KEYNOTES

A1 EXTENT OF EXISTING DEEPER STRUCTURAL SLAB. EXACT LOCATION TO BE CONFIRMED IN FIELD, SEE STRUCTURAL DRAWINGS FOR ADDITIONA INFORMATION.

A3 EXISTING FIRE EXTINGUISHER CABINET TO REMAIN. CONFIRM MAKE AND MODEL IS MRI SAFE.

A4 PROVIDE KEYLOCK ON EXISTING DOOR. A5 DOOR BY MRI SHEILDING VENDOR

BY RECESSED COUNTERTOP BRACKETS.

A2 PROVIDE MRI SAFE ADJUSTABLE BLINDS.

A6 3FORM ILTUO FLOOR TO CEILING 200.58 PARTITION. REFER TO FINISH PLAN FOR RESIN FINISH. FASTEN TO FLOOR AND CEILING, PROVIDE BLOCKING ABOVE CEILING AT VERTICAL SUPPORTS.

A7 NEW SHAFT FOR QUENCH VENT. MORE INFORMATION PROVIDED ON MEP DRAWINGS AND OPENING INFORMATION ON STRUCTURAL DRAWINGS. COORDINATE WORK SCHEDULE WITH CURRENT OCCUPANTS OF SPACE. A8 ADD ALTERNATE: DEMOLISH EXISTING PLASTIC LAMINATE COUNTERTOP. REPLACE WITH NEW SOLID SURFACE COUNTERTOP, SSP-1, SUPPORTED

A9 PROVIDE MINIMUM 1" RADIUS ON SOLID SURFACE EDGE A10 PROVIDE 24" X 60" STAINLESS STEEL FRAMED WALL MOUNTED MIRROR. PROVIDE IN-WALL BLOCKING.

A11 APPROXIMATE CONCRETE INFILL OF TRENCH FOR SANITARY LINE FOR FLOOR DRAIN. A12 CONFIRM EXACT LOCATION OF EXISTING OUTLET, MAY NEED TO BE RELOCATED TO ACCOMMODATE MILLWORK. REVIEW WITH ENGINEER

A13 RELOCATED FEC.

A14 LOCKABLE UPPER CABINET A15 PATCH WALL AND FLOOR AS TO LEAVE NO EVIDENCE OF REMOVAL. MATCH EXISTING FINISHES. ADD CEILING TILES AND GRID IN AREA OF WALL REMOVAL TO MATCH EXISTING.

A16 GC TO PROVIDE AND INSTALL MRI IN-USE LIGHT. PROVIDE DOOR CONTACT TO COORDINATE MRI USE. A17 SLAB PENETRATIONS TO BE 2HR RATED

A18 REMOVABLE ACCESS RF WALL FOR MRI RIGGING A19 GROMMET OPENING FOR INJECTOR CABLING. COORDINATE WITH **VENDOR DOCUMENTS** A20 EPO BUTTON. COORDINATE WITH VENDOR.

A21 PROVIDE 2.5" DIA. WAVEGUIDE PIPE FOR INJECTOR CONTROL. PROVIDE COVER FOR WHEN THE WAVEGUIDE PIPE IS NOT IN USE. A22 RECESSED COUNTER SUPPORT BRACKETS. BASIS OF DESIGN: RAKKS

E-1824 INSIDE WALL MOUNT A23 ALL OUTLETS PROVIDED IN NEW MRI ROOM - W189C TO BE MRI SAFE/ RF

A24 PROVIDE MONITOR BRACKETS A25 GC TO INSTALL HOLLOW METAL DOOR AND FRAME TO KEEP MRI 1 OPERATIONAL DURING THE COURSE OF CONSTRUCTION. DOOR TO BE A

RIGGING PLAN.

SMOKE DOOR WITH GASKETING. TOWARD THE END OF CONSTRUCTION, GC TO REMOVE DOOR AND FRAM, RESTORE WALL TO ORIGINAL CONDITION AND RATING.

A26 GC TO REINSTALL SYSTEMS FURNITURE AFTER RIGGING. REFER TO

STRUCTURAL ENGINEER

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

212-689-3110

CONSULTANTS:

MEP ENGINEER

LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001 PHONE: 646.674.6100

ARRAY-ARCHITECTS.COM

NEW YORK, NY. 10016

470 PARK AVE SOUTH, 11th FLOOR

Architect of Record:

OWNER:

COLUMBIA DOCTOR'S TARRYTOWN

PROJECT: **NEW MRI**

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591

KEY PLAN

W - WEST BUILDING **E - EAST BUILDING** NORTH NORTH

DESCRIPTION DATE **REVISIONS/ISSUES**

SHEET TITLE: INTERIOR ELEVATIONS

SEAL:

CON/REF No. CONTRACT No. **SCALE:** AS NOTED PROJECT No. 6109 CHECKED: CH DRAWN: KU


DATE:06/18/2021

A-600.00

DWG OF 47

SHEET NO.

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be prosecuted under the law.



PARTITION TYPES GENERAL NOTES

Architect of Record:

ARRAY-ARCHITECTS.COM 470 PARK AVE SOUTH, 11th FLOOR

CONSULTANTS:

212-689-3110

NEW YORK, NY. 10016

STRUCTURAL ENGINEER

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER

PHONE: 646.674.6100

LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001

OWNER:

COLUMBIA DOCTOR'S **TARRYTOWN**

PROJECT: **NEW MRI**

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591

KEY PLAN

W - WEST BUILDING **E - EAST BUILDING**

EXISTING PARTITION TYPE LEGEND

EXISTING SMOKE PARTITION TO REMAIN. CONTRACTOR TO INSPECT AND SURVEY EXISTING CONDITION OF PARTITION AND PERFORM ALL NECESSARY WORK REQUIRED TO CONFORM WITH PARTITION TYPE *. CONTRACTOR TO USE MATERIALS IDENTICAL TO EXISTING MATERIALS OR MATERIALS WHOSE INSTALLED PERFORMANCE WILL EQUAL OR SURPASS THAT OF PARTITION TYPE * MATERIALS.

EXIST. SMOKE PARTITION

EXISTING PARTITION TO REMAIN WITH LEAD LINING ON SIDE OF PARTITION WHERE RADIATION EMITTING EQUIPMENT IS LOCATED. CONTRACTOR TO INSPECT AND SURVEY EXISTING CONDITION OF PARTITION AND PERFORM ALL NECESSARY WORK TO ACHIEVE THE REQUIREMENTS SET FORTH IN THE PHYSICIST REPORT. RADIATION TEST OF EXISTING ROOM IS REQUIRED TO DETERMINE ITS COMPLIANCE. IF REQUIREMENTS ARE NOT MET, CONTRACTOR TO MODIFY PARTITION FROM CEILING TO DECK TO CONFORM WITH PARTITION TYPE * CONTRACTOR TO USE MATERIALS IDENTICAL TO EXISTING MATERIALS OR MATERIALS WHOSE INSTALLED PERFORMANCE WILL EQUAL OR SURPASS THAT OF PARTITION TYPE * MATERIALS.

EXIST. LEAD LINED PARTITION

CONTRACTOR TO MODIFY EXISTING PARTITION FROM CEILING TO DECK TO CONFORM WITH PARTITION TYPE */UL #465. CONTRACTOR, ALSO, TO INSPECT AND SURVEY EXISTING CONDITION OF PARTITION AND PERFORM ALL WORK TO MAINTAIN REQUIRED FIRE RATING FROM FLOOR TO STRUCTURE ABOVE. CONTRACTOR TO USE MATERIALS IDENTICAL TO EXISTING MATERIALS OR MATERIALS WHOSE INSTALLED PERFORMANCE WILL EQUAL OR SURPASS THAT OF PARTITION TYPE */UL #465 MATERIALS.

EXIST. PARTITION MODIFIED TO A 2 HOUR

DESCRIPTION DATE **REVISIONS/ISSUES**

SHEET TITLE: PARTITION TYPES & RATED ASSEMBLY

SEAL:

DATE:06/18/2021 CON/REF No. CONTRACT No. **SCALE:** AS NOTED **PROJECT No.** 6109 CHECKED: CH DRAWN: KU


A-720.00

C COPYRIGHT - 2020 ARRAY-ARCHITECTS P.C.

DWG OF 47

SHEET NO.

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be prosecuted under the law.



Architect of Record: ARRAY-ARCHITECTS.COM 470 PARK AVE SOUTH, 11th FLOOR **NEW YORK, NY. 10016**

CONSULTANTS:

STRUCTURAL ENGINEER

PHONE: 570.496.7020

MEP ENGINEER

LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001

OWNER:

COLUMBIA DOCTOR'S TARRYTOWN

PROJECT: NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591

DESCRIPTION DATE **REVISIONS/ISSUES**

SHEET TITLE: FRAMING DETAILS

SEAL:

CON/REF No. CONTRACT No. SCALE: AS NOTED PROJECT No. 6109 CHECKED: Checker **DRAWN:** Author

DATE:06/18/21


SHEET NO. A-721.00


C COPYRIGHT - 2020 ARRAY-ARCHITECTS P.C.


DWG OF 47


This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be

prosecuted under the law.









F2 PROTECT EXISTING FINISHES TO REMAIN DURING CONSTRUCTION.

MATCH EXISTING PAINT OR INDICATED PAINT ACCENT PAINT COLOR.

COORDINATE WITH FACILITY AND EQUIPMENT TO PERFORM WORK.


F3 FLASH PATCH FEATHERED UP TO MEET EDGE OF RF SHIELDING.

EXISTING MRI WILL HAVE OT BE RAMPED DOWN TO INSTALL.

F6 PATCH AND MATCH EXISTING FLOORING WHERE REQUIRED.

PATCH AND REPAIR AS NECESSARY.

F5 INSTALL SALVAGED CORNER GUARD.



A6 FIRST FLOOR FINISH AND FLOOR PATTERN PLAN
1/4" = 1'-0"

ARRAY-ARCHITECTS.COM 470 PARK AVE SOUTH, 11th FLOOR **NEW YORK, NY. 10016** 212-689-3110

Architect of Record:

CONSULTANTS:

STRUCTURAL ENGINEER

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER


LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001 PHONE: 646.674.6100

OWNER:

COLUMBIA DOCTOR'S TARRYTOWN

PROJECT: NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591



DESCRIPTION REVISIONS/ISSUES

SHEET TITLE: FINISH PLAN AND SCHEDULE

SEAL:

DATE:06/18/2021 CON/REF No. CONTRACT No. SCALE: AS NOTED PROJECT No. 6109 CHECKED: CH DRAWN: KU


SHEET NO.

A-800.00

C COPYRIGHT - 2020 ARRAY-ARCHITECTS P.C.

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be prosecuted under the law.

CD Submission 06/18/2021



ARRAY-ARCHITECTS.COM 470 PARK AVE SOUTH, 11th FLOOR NEW YORK, NY. 10016

Architect of Record:

CONSULTANTS:

STRUCTURAL ENGINEER

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER

LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001

COLUMBIA DOCTOR'S TARRYTOWN

NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591

W - WEST BUILDING **E - EAST BUILDING**

DESCRIPTION REVISIONS/ISSUES

SHEET TITLE: DOOR TYPES, FRAME DETAILS, SCHEDULE & BORROWED LIGHT DETAILS, SCHEDULE

DATE:06/18/2021 CON/REF No. CONTRACT No. **SCALE:** AS NOTED PROJECT No. 6109 CHECKED: CH DRAWN: KU

A-900.00

COPYRIGHT - 2020 ARRAY-ARCHITECTS P.C.

DWG OF 47

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be

CONSULTANTS:

STRUCTURAL ENGINEER

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER

PHONE: 646.674.6100


LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001


OWNER:

COLUMBIA DOCTOR'S TARRYTOWN

PROJECT: NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591





DESCRIPTION DATE **REVISIONS/ISSUES**

SHEET TITLE: GE MRI-C1

SEAL:


CON/REF No. CONTRACT No. SCALE: AS NOTED PROJECT No. 6109 CHECKED: CH DRAWN: KU

DATE:06/18/2021

SHEET NO. V-101.00

DWG OF 47

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be prosecuted under the law.



DISCLAIMER

GENERAL SPECIFICATIONS

- GE is not responsible for the installation of developers and associated equipment, lighting, cassette trays and protective screens or derivatives not mentioned in the order.
- The final study contains recommendations for the location of GE equipment and associated devices, electrical wiring and room arrangements. When preparing the study, every effort has been made to consider every aspect of the actual equipment expected to be installed.
- . The layout of the equipment offered by GE, the dimensions given for the premises, the details provided for the pre-installation work and electrical power supply are given according to the information noted during on-site study and the wishes expressed by the customer.
- The room dimensions used to create the equipment layout may originate from a previous layout and may not be accurate as they may not have been verified on site. GE cannot take any responsibility for errors due to lack of information.
- · Dimensions apply to finished surfaces of the room.

Columbia University Medical Center

- Actual configuration may differ from options presented in some typical views or tables.
- · If this set of final drawings has been approved by the customer, any subsequent modification of the site must be subject to further investigation by GE about the feasibility of installing the equipment. Any reservations must be noted.
- The equipment layout indicates the placement and interconnection of the indicated equipment components. There may be local requirements that could impact the placement of these components. It remains the customer's responsibility to ensure that the site and final equipment placement complies with all applicable
- · All work required to install GE equipment must be carried out in compliance with the building regulations and the safety standards of legal force in the country concerned.
- These drawings are not to be used for actual construction purposes. The company cannot take responsibility for any damage resulting therefrom.

CUSTOMER RESPONSIBILITIES

- It is the responsibility of the customer to prepare the site in accordance with the specifications stated in the final study. A detailed site readiness checklist is provided by GE. It is the responsibility of the customer to ensure all requirements are fulfilled and that the site conforms to all specifications defined in the checklist and final study. The GE Project Manager of Installation (PMI) will work in cooperation with the customer to follow up and ensure that actions in the checklist are complete, and if necessary, will aid in the rescheduling of the delivery and installation date.
- Prior to installation, a structrual engineer of record must ensure that the floor and ceiling is designed in such a way that the loads of the installed system can be securely borne and transferred. The layout of additional structural elements, dimensioning and the selection of appropriate installation methods are the sole responsibility of the structural engineer. Execution of load bearing structures supporting equipment on the ceiling, floor or walls are the customer's responsibility.

	TIFIES THAT I HAVE READ AND APPROVED	
DATE	NAME	SIGNATURE

SIGNA ARTIST

GLOBAL SITE READINESS CHECKLIST (DI)

DOC1809666 Rev. 7

Site Ready Checks at Installation

EHS Site Requirements

Enough space to store tools, equipment, parts, install waste and the general area free from obstruction and trip hazards.

Enough necessary facilities for the GE employees available.

No 3rd parties working in the area that may affect the safety of the installation activity.

Area free from any chemical, gas, dust, welding fume exposure and has painting been completed and dry.

All emergency routes identified, signed and clear from obstruction.

Overall access route to the scan room free from obstruction / high hazards

Accessible single source lockable panel that LOTO can be applied to for GE equipment installation (MDP and/or PDU).

There are no other conditions or hazards that you have observed or have been made aware of by the customer or contractors on site.

Specific for MR Magnet Delivery

Site Ready Checks for MR Magnet delivery if separate from system

All permits, plans and permissions received for rigging and/or delivery.

Adequate delivery route from truck to final place of installation has been reviewed with all stakeholders, all communications/notifications have occurred, arrangements have been made for special handling (rigging, elevator, fork lift, etc.)

All floors along delivery route will support weight of the equipment, temporary reinforcements arranged if needed.

Chilled water supply for Water Cooled Compressor or Air Cooled Compressor is ready and meets GE specifications.

Water drain available in the equipment room, if applicable.

Power for MR compressor & Chiller is available.

Power is available for magnet monitoring.

Connectivity is available for magnet monitoring. Ensure cryogen venting system is complete to outside the building & available for magnet connection.

Exhaust fan system is installed per GE requirements and VERIFIED operational by suction at intake.

Required for Mechanical Install start

Room dimensions, including ceiling height, for all Exam, Equipment/Technical & Control rooms meets GE specifications. Ceiling support structure, if indicated on the GE drawing, is in the correct location and at the correct height according to the Original

Equipment Manufacturer specifications.

Levelness and spacing has been measured, and is ready for the installation of any GE supplied components.

Finished ceiling is installed. If applicable ceiling tiles installed per PMI discretion.

Floor levelness/flatness is measured and within tolerance, and there are no visible defects per GEHC specifications.

Entry door threshold meets PIM requirement

Floor Strength and thickness have been discussed with customer/contractor and they have confirmed GE requirements are met.

Rooms that will contain equipment, including staging areas if applicable, are construction debris free. Precautions must be taken to prevent debris from entering rooms containing equipment.

Cable ways (floor/wall/ceiling/Access Flooring) are available for installation of GE cables are of correct length and diameter.

Overhead cable trays must be install exactly as shown in GE supplied installation draings to assure proper Gradient Cable length.

Cable ways routes per GE Final drawings and cable access openings areas installed at a time determined by GEHC PM. Surface floor duct can be installed at time of system installation.

Adequate room illumination installed and working.

Customer supplied countertops where GE equipment will be installed are in place.

RF Shield installed with possible exception of magnet entrance. RF Shield Effectivity and Ground Isolation Test needed. If GE responsible for supplying RF shield, the RF shield Effectivity and Ground Isolation Test data is a Mandatory attachment into MyProjects.

Required for Calibration start

HVAC systems Installed, and the site meets minimum environmental operational system requirements.

System power & grounding (PDB/MDP) is available as per GE specifications.

System power & grounding (PDB/MDP) is installed at point of final connection and ready to use. Lock Out Tag Out is available.

PMI to confirm all feeder wires and breaker are size appropriately. EPO installed if needed. PMI to confirm with electrician all power and signal cables are well terminated ensuring there are no loose connections.

Network outlets installed.

MRI-M221448-FIN-00.DWG

Computer network available and working.

Note: The details shown here are only an extract from DOC1809666. For the complete document please contact your PMI.

- |Rev A|Date 30/JUL/2020 |

C1 - Cover Sheet

02/27

MEP ENGINEER

LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001 PHONE: 646.674.6100

ARRAY-ARCHITECTS.COM

STRUCTURAL ENGINEER

NEW YORK, NY. 10016

212-689-3110

CONSULTANTS:

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

470 PARK AVE SOUTH, 11th FLOOR


Architect of Record:


OWNER:

COLUMBIA DOCTOR'S **TARRYTOWN**

PROJECT: **NEW MRI**

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591





DESCRIPTION REVISIONS/ISSUES

SHEET TITLE: GE MR-C2

SEAL:

DATE:06/18/2021 CON/REF No. CONTRACT No. **SCALE:** AS NOTED **PROJECT No.** 6109 CHECKED: CH DRAWN: KU

SHEET NO. V-102.00

DWG OF 47

CD Submission 06/18/2021

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be prosecuted under the law.

470 PARK AVE SOUTH, 11th FLOOR

ARRAY-ARCHITECTS.COM

NEW YORK, NY. 10016

CUSTOMER SITE READINESS REQUIREMENTS

- Any deviation from these drawings must be communicated in writing to and reviewed by your local GE healthcare installation project manager prior to making changes.
- Make arrangements for any rigging, special handling, or facility modifications that must be made to deliver the equipment to the installation site. If desired, your local GE healthcare installation project manager can supply a reference list of rigging contractors.
- New construction requires the following;
- Secure area for equipment,
- Power for drills and other test equipment,

Capability for image analysis,

- Restrooms.
- Provide for refuse removal and disposal (e.g. crates, cartons, packing)
- It is the customer's responsibility to contract a vibration consultant/engineer to implement site design modifications to meet the GE vibration specification. Refer to the system preinstallation manual for the vibration specification.

MRI SITE PLANNING REMINDERS

Please refer to pre-installation checklist in pre-installation manual listed on the cover sheet for items critical to image quality.

- The layout should be arranged so that the 5g line is contained to the magnet room. If not possible, a barrier is recommended to prevent entry to the 5g field area.
- 2. The spaces around, above, and below the magnet must be reviewed for effects of the 5g, 3g, 1g, and .5g fields. Refer to the proximity limit chart in the MR pre-installation manual referenced on the cover sheet.
- 3. For moving metal, the restriction lines typically extend outside of the MR space. Please confirm there are no moving metal concerns within these areas.
- 4. For vibration, analysis to be completed as required per pre-installation manual.
- 5. For EMI, review the site for the location of the main electrical feeders, AC devices, or distribution systems. An EMI study is recommended if large AC systems are nearby.
- 6. Details of the floor below the magnet must be reviewed. The structural engineer must verify that the quantity of steel in the volume 10ft [3.1m] x 10ft [3.1m] x 1ft [.3m] deep (below the magnet) does not exceed the allowable steel content as given in the MR pre-installation manual referenced on the cover sheet.
- Remove, cover, or fill-in abandoned ducts or troughs from the Equipment and Magnet rooms. Access/computer room flooring in the Equipment room can either be removed or assessed and reinforced to support heavier cabinets.

Responsibility for the coordination, design, engineering, and site preparation resides with the customer and their project architects and contractors. GE does not, by providing reviews and furnishing comments and assistance, accept any responsibility beyond its obligations as defined in the MR system, sale/purchase agreement.

IMAGE QUALITY CONSIDERATIONS

Broadband RF noise is a single transient or continuous series of transient disturbances caused by an electrical discharge. Low humidity environmental conditions will have higher probability of electrical discharge. The electrical discharge can occur due to electrical arcing (micro arcing) or merely static discharge. Some potential sources capable of producing electrical discharge include:

- Loose hardware/fasteners vibration or movement (electrical contunuity must always be maintained)
- Flooring material including raised access flooring (panels & support hardware) and carpeting
- Electrical fixtures (i.e. Lighting fixtures, track lighting, emergency lighting, battery chargers, outlets)
- Ducting for HVAC and cable routing
- RF shield seals (walls, doors, windows etc.)

For additional information regarding image quality, refer to the pre-installation manual listed on the cover sheet.

MAGNETIC INTERFERENCE SPECIFICATIONS

- The customer must establish protocols to prevent persons with cardiac pacemakers, neurostimulators, and biostimulation devices from entering magnetic fields of greater than 5 gauss (exclustion zone).
- Main power transformers must remain outside the 3 gauss field. EMI < 20mg rms ac. EMI < 5.87mg dc. Potential exists under fault conditions that the 5 gauss line may expand radially to 9.35 ft. [2.85 m] and axially
- to 14.27 ft. [4.35 m] for 1 seconds or less. It should be noted that normal rampdowns or magnet rundown unit initiated quenches will not cause the magnetic field to expand.
- It is recommended every site consider the event of a quench and plan accordingly (such as placing 5 gauss warning signs at expanded locations).
- The ferrous metal objects listed below must not move into or inside of the moving metal sensitivity line during scans.

TYPCIAL MOVING MAGNETIC MASS	DISTANCE	RADIALLY	DISTANCE AXIALLY		
Carts, Gurneys 100-400 lbs [45-182 kg]	3 Gau	ss line	3 Gauss line		
Forklifts, small elevator, cars, minivans vans, pickup trucks, ambulances (objects greater than 400 lbs [182 kg])	15.5 FT	4.72 M	24.6 FT	7.5 M	
Buses and trucks (dump, tractor trailer, utility, fire trucks)	18.1 FT	5.52 M	28.75 FT	8.76 N	


03/27 Columbia University Medical Center MRI-M221448-FIN-00.DWG Rev AlDate 30/JUL/2020 SIGNA ARTIST C1 - Cover Sheet


OWNER:

COLUMBIA DOCTOR'S TARRYTOWN

PROJECT: NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591





DESCRIPTION REVISIONS/ISSUES

SHEET TITLE: GE MR-A1

SEAL:

DATE:06/18/2021 CON/REF No. CONTRACT No. **SCALE:** AS NOTED **PROJECT No.** 6109 CHECKED: CH DRAWN: KU

SHEET NO. V-103.00

DWG OF 47

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be prosecuted under the law.

MAX HEAT MAX HEAT MAX HEAT LEGEND WEIGHT OUTPUT OUTPUT OUTPUT DESCRIPTION DESCRIPTION BY ITEM DESCRIPTION A GE Supplied D Available from GE A 14 Operator workspace GE Supplied/contractor A 1 1.5T Magnet 8189 26 Minimum opening for equipment delivery is 40 in. w x 82 in. Existing h, contingent on a 72 in. corridor width 2 Rear pedestal A 15 | Pneumatic patient alert 0.5 Customer/contractor ---- 5 Gauss A 16 Operators console desk 3 GEM Patient table 125 supplied and installed Minimum opening for equipment delivery is 43 in. w x 82 in. h, contingent on a 96 in. corridor width D 17 MR Elastography 3, 1, 0.5 Gauss 200 Gauss 4 Magnet rundown unit 480 53.4 Magnetic Shielding 100, 50, 30, 10 Gauss A 18 Music system 5 Phantom set storage cabinet Counter top for equipment- provide grommeted openings as required to route cables 6 Blower box 1535 A 19 Injector on pedestal 94 RF SHIELD - 100 dB ATTENUATION A 20 Injector control 675 17 Base cabinet for storage of: surface coils, patient positioning 7 Penetration cabinet 10697 pads, phantoms, etc. A 21 Injector power supply 660

130

4000

901


C 31 Louvered doors - refer to preinstall for requirements

B 22 Main disconnect panel

B 23 Dimplex chiller

A 8 Secondary penetration wall

A 9 Power, gradient, RF cabinet



Architect of Record: architects

ARRAY-ARCHITECTS.COM 470 PARK AVE SOUTH, 11th FLOOR **NEW YORK, NY. 10016** 212-689-3110

CONSULTANTS:

STRUCTURAL ENGINEER

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER

LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001 PHONE: 646.674.6100

OWNER:

EXAM ROOM HEIGHT


13'-0"

FINISHED FLOOR TO SLAB HEIGHT

COLUMBIA DOCTOR'S TARRYTOWN

PROJECT: NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591



W - WEST BUILDING E - EAST BUILDING KEY PLAN

DESCRIPTION REVISIONS/ISSUES

SHEET TITLE: GE MR-A2

SEAL:

DATE:06/18/2021 CON/REF No. CONTRACT No. SCALE: AS NOTED PROJECT No. 6109 CHECKED: CH DRAWN: KU

SHEET NO. V-104.00

DWG OF 47

Note: Define RF shield's inset according to provisions made by the RF Shield vendor.

CONTROL ROOM

MAGNET ROOM

ARRAY-ARCHITECTS.COM 470 PARK AVE SOUTH, 11th FLOOR NEW YORK, NY. 10016 212-689-3110

CONSULTANTS:

STRUCTURAL ENGINEER

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER

SIDE VIEW WITH MAGNETIC FIELD

Cable tray

-Unistrut


LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001 PHONE: 646.674.6100


OWNER:

COLUMBIA DOCTOR'S TARRYTOWN

PROJECT:
NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591







SHEET TITLE:
GE MR-A3

SEAL:

DATE: 06/18/2021

CON/REF No.

CONTRACT No.

SCALE: AS NOTED

CON/REF No.
CONTRACT No.
SCALE: AS NOTED
PROJECT No. 6109
CHECKED: CH
DRAWN: KU


V-105.00

DWG OF 47

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be

C COPYRIGHT - 2020 ARRAY-ARCHITECTS P.C.

CD Submission 06/18/2021
This document is a copyrigh property of Array Architects project only. Reproduction of written permission of Array prosecuted under the law.



Ε

_

С

2:38:14 PM 12" = 1'-0"

N K magnet's isocenter.

Architect of Record:

ARRAY-ARCHITECTS.COM 470 PARK AVE SOUTH, 11th FLOOR **NEW YORK, NY. 10016** 212-689-3110

CONSULTANTS:

STRUCTURAL ENGINEER

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER


LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001 PHONE: 646.674.6100


OWNER:

COLUMBIA DOCTOR'S TARRYTOWN

PROJECT: NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591







SHEET TITLE: GE MRI-A4

SEAL:

CON/REF No. CONTRACT No. SCALE: AS NOTED PROJECT No. 6109 CHECKED: CH DRAWN: KU

DATE:06/18/2021

SHEET NO. V-106.00

DWG OF 47

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be prosecuted under the law.

C COPYRIGHT - 2020 ARRAY-ARCHITECTS P.C.


ACOUSTICS SPECIFICATIONS


Acoustic and vibroacoustic information is provided for site planning and architectural design activities. It is the customer's responsibility to hire a qualified acoustic engineer for solutions to further attenuate this transmitted noise and vibration, if required. The actual room noise level may vary based on room design, optional equipment, and usage:

0.4

Control Room: 62 dBA Equipment Room: 80 dBA Magnet Room: 127 dBA* (maximum sound pressure level at magnet bore isocenter)

* Frequency: 20 Hz to 20kHz





FREQUENCY [Hz]

SIGNA ARTIST

Low Frequency Magnet Floor Vibration (Vibration Amplitude at Each Foot)

MAGNETIC PROXIMITY LIMITS

ISOGAUSS PLOTS

The isogauss contour plots illustrated represent the performance of the system and the steel room shield

described in the final shield design dated JUL. 27, 2020. The actual magnetic flux density at any point in the

contours shown are only approximations of actual flux densities found at corresponding distances from the

vicinity of the magnet when installed may vary from the contour plots due to factors such as the concentrating

effects of nearby ferrous objects and ambient magnetic fields, including the earth's magnetic field. Therefore, the


Gauss (mT) Limit	Equipment
0.5 gauss (0.05mT)	Nuclear camera
1 gauss (0.1mT)	Positron Emission Tomography scanner, Linear Accelerator, Cyclotrons, Accurate measuring scale, Image intensifiers, Bone Densitometers, Video display (tube), CT scanner, Ultrasound, Lithotriptor, Electron microscope, Digital X-Ray
3 gauss (0.3mT)	Power transformers, Main electrical distribution transformers
5 gauss (0.5mT)	Cardiac pacemakers, Neurostimulators, Biostimulation devices
10 gauss (1mT)	Magnetic computer media, Line printers, Film processor, X-ray tubes, Emergency generators, Commercial laundry equipment, Food preparation area, Water cooling equipment, HVAC equipment, Major mechanical equipment room, Credit cards, watches, and clocks, Air conditioning equipment, Fuel storage tanks, Motors greater than 5 horsepower
50 gauss (5mT)	Metal detector for screening, LCD panels, Telephones
No Limit	Digital Detectors

The customer must provide detail defining ferrous material below the magnet to the Project Manager so the GE Healthcare MR Siting and Shielding team can review for compliance.

Limits Of	Steel Mass	Distance From Ma	ignet Isocenter	Distance Below Top	Surface Of Floor
kg/m²	lbs/ft²	mm	in	mm	in
0	0	0 - 1143	0-45	0 - 76	0-3
9.8	2	1143 - 1194	45-47	76 - 127	3-5
14.7	3	1194 - 1321	47-52	127 - 254	5-10
39.2	8	1321 - 1397	52-55	254 - 330	10-13
98.0	20	1397+	55+	330+	13+

The actual field strength can be affected by Magnetic shielding, Earth's magnetic field, other magnetic fields and stationary or moving metal. This information must be used to evaluate potential site interaction of GE Healthcare equipment with other non-GE Healthcare equipment. Magnetic shielding can be installed to prevent interaction between the magnet and nearby sensitive devices. The GE Healthcare Project Manager of Installation (PMI) can work with the customer to coordinate the magnetic shielding site evaluation. The customer is responsible for installation of all magnetic shielding. Rev AlDate 30/JUL/2020 C1 - Cover Sheet 06/27 MRI-M221448-FIN-00.DWG

SOUND PRESSURE SPECTRAL DISTRIBUTION



Columbia University Medical Center

MRI-M221448-FIN-00.DWG

| 1:10 |Rev A|Date 30/JUL/2020

ARRAY-ARCHITECTS.COM 470 PARK AVE SOUTH, 11th FLOOR **NEW YORK, NY. 10016** 212-689-3110

CONSULTANTS:

STRUCTURAL ENGINEER

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER

LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001 PHONE: 646.674.6100


OWNER:

COLUMBIA DOCTOR'S TARRYTOWN

PROJECT: NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591





DATE **DESCRIPTION REVISIONS/ISSUES**

SHEET TITLE: GE MR-A5

SEAL:

CON/REF No. CONTRACT No. **SCALE:** AS NOTED PROJECT No. 6109 CHECKED: CH DRAWN: KU

DATE:06/18/2021


SHEET NO. V-107.00

DWG OF 47

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be prosecuted under the law.

C COPYRIGHT - 2020 ARRAY-ARCHITECTS P.C.





SIGNA ARTIST

07/27

RF shielding

SCALE 1:30

Columbia University Medical Center

CONSULTANTS:

STRUCTURAL ENGINEER

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER


LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001 PHONE: 646.674.6100

OWNER:

COLUMBIA DOCTOR'S TARRYTOWN

PROJECT: **NEW MRI**

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591





SHEET TITLE: GE MRI-A6

KEY PLAN


SEAL: **DATE:**06/18/2021

CON/REF No. CONTRACT No. SCALE: AS NOTED PROJECT No. 6109 CHECKED: CH DRAWN: KU

SHEET NO. V-108.00

DWG OF 47

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be prosecuted under the law.



CONSULTANTS:

STRUCTURAL ENGINEER

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER


LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001 PHONE: 646.674.6100


OWNER:

COLUMBIA DOCTOR'S TARRYTOWN

PROJECT: **NEW MRI**

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591







SHEET TITLE: GE MRI-A7

SEAL: **DATE:**06/18/2021 CON/REF No.


CONTRACT No. SCALE: AS NOTED PROJECT No. 6109 CHECKED: CH DRAWN: KU

SHEET NO. V-109.00

C COPYRIGHT - 2020 ARRAY-ARCHITECTS P.C.

DWG OF 47

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be prosecuted under the law.





CONSULTANTS:

STRUCTURAL ENGINEER

REUTHER + BOWEN DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER

PHONE: 646.674.6100

LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001

CRITICAL ITEMS FOR MAGNET DELIVERY

ROUTING


- The customer is solely liable for routing of components from dock to final site.
- GE must be able to move system components in or out with no need to uncrate or disassemble any of the components. The entire passageway must be cleared, adequately lighted and free from dust.

DELIVERY

- The floor and it surfacing must be able to withstand the live load of components and handling equipment.
- Floor surfacing must be continuous.
- The customer must protect any fragile flooring surfaces.

MINIMUM SPECIFICATIONS FOR MAGNET ROUTING

- Floor must be able to withstand a moving load of 4823kg
- Height: 2238 mm [88.11in], width: 2303 mm [90.66in]



INSTALLATION AND DELIVERY ACCEPTANCE

If the site is not ready, GE can delay the delivery time.

Columbia University Medical Center

decision for the delivery time.

A survey of the site established by the customer and GE will make the

. This survey of the site (a form is made available by GE) is only to check if


the apparent conditions of the site allow the equipment to be delivered.

2303 [90.7 in]

MRI-M221448-FIN-00.DWG

SIGNA ARTIST


SIDE VIEW OF MAGNET FRONT VIEW OF MAGNET Recommended opening for side (wall) delivery: 2400 mm [90.55in] (width) x 2500 mm [98.42in] (height)



24/7 chilled water and 480v power for shield/cryo cooler

- 24/7 120v power for the magnet monitor
- Phone lines for magnet monitoring and emergency use
- Magnet room exhaust fan
- Cryogen venting (if roof hatch, completed within 24 hrs)

This is only a partial list of items required for delivery of the magnet. For a complete checklist refer to the pre-installation manual referenced on cover sheet.



Front view of the BRM Gradient

DIMENSIONS

Side view of the BRM Gradient

C1 - Cover Sheet

EQUIPMENT	LxW	xH	WE	IGHT	NOTE
B-Sessione Val	mm	in	kg	lbs	
Replacement BRM gradient coil assembly on a shipping cradle/cart	991x2536x1499	39x99.84x59	1449	3194	Initial gradient coil assembly is shipped installed in the magnet. Shipping/installation cart is used to install re-placement coil assembly only.

The weight bearing structure of the site should support any additional weight of the main replacement parts occurring during maintenance of the magnet, throughout the whole lifecycle of the MR.


Rev AlDate 30/JUL/2020


OWNER:

TARRYTOWN PROJECT: NEW MRI

COLUMBIA DOCTOR'S

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591





DESCRIPTION REVISIONS/ISSUES

SHEET TITLE: GE MRI-A8

SEAL:

DATE:06/18/2021 CON/REF No. CONTRACT No. **SCALE:** AS NOTED PROJECT No. 6109 CHECKED: CH DRAWN: KU

SHEET NO. V-110.00


DWG OF 47

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title

10/27

CD Submission 06/18/2021

SHIELDING PLAN [REFERENCE ONLY-NOT TO SCALE]



All surfaces are 1/8" M36

ADDITIONAL NOTES

- Attractive force of a magnet on the ceiling and walls of the steel shield will be a maximum of 30% of the shield's weight. Support structure must account for these additional forces.
- Attractive forces on the floor or sub-floor may be sufficient to lift shield's weight. Support structure must prevent this with adhesive or bracing.
- M36 Silicon steel is used in this shield. see shielding requirements sheets for specifications.

Columbia University Medical Center

PERFORMANCE STATEMENT FOR MAGNETIC SHIELD DESIGN

GE HEALTHCARE ("GE") provides this magnetic shielding design solely for use with GE's Magnetic Resonance System ("System"). This

The magnetic shield, shown in this set, when constructed in strict accordance with this magnetic shield design will:

("Design Objectives"), subject to and conditioned upon the following:

- The magnetic field strength measurements are conducted using techniques which account for the effects of the earth's
- 2. The magnetic field strength measurements are conducted using techniques which account for magnetic field distortions caused by magnetic materials in the vicinity of the measurement and;
- specifications and:
- the System has not been and is not being subjected to improper maintenance, and (iii) Customer has complied with all applicable instructions or recommendations of GE related to the System and magnetic shield; and
- Customer operates the System at or below its specified field strength.

GE will re-perform the magnetic shielding design services, to the extent that GE determines the design services fail to meet the Design Objectives specified above, provided GE is notified in writing by Customer of such failure within one (1) week after the date of availability of the System for first use, subject to all other terms and conditions of the agreement between GE and Customer for the purchase or lease of the System.

The foregoing sets forth Customer's exclusive remedies and GE's sole liability for claims based on the failure of the magnetic shielding design services to meet the Design Objectives. In no event shall GE be liable for special or consequential damages.

The magnetic field gradients resulting from the magnetic shield's effect on the System's magnet will allow the System's magnet to meet published homogeneity specifications upon completion of required installation procedures.

The final performance of the magnetic shield set forth in this drawing set depends on strict adherence to the construction and material requirements detailed herein. Any deviation from these requirements will degrade shielding effectiveness. For this reason, GE recommends that only established magnetic shielding vendors be used. Vendors should be required to provide documentary evidence to demonstrate that their shielding product complies with GE's construction requirements, material requirements, chemical composition specifications, and annealing specifications; or the material's DC magnetization curve at high induction is equivalent to GE's specification on SH2.. A guarantee that the finished shield will meet the Customer's fringe field containment requirements should be sought from the shielding vendor.

The isogauss contours shown on this drawing set indicate the magnetic fringe field containment which corresponds to the predicted interaction between the System and the magnetic shield design described. GE creates these isogauss contour plots only to predict if the shield design meets the Design Objectives. All other effects such as the superposition of the earth's magnetic field, residual magnetism, and localized field concentration effects, due to ferromagnetic structural elements, are not considered or indicated. The exact location of a particular isogauss contour may differ from that shown due to: (1.) the above stated reasons and (2.) the resolution limits of the mathematical modeling techniques used to derive these results. Note also that for actively shielded Systems, fault conditions may exist where the isogauss contours expand temporarily (refer to your System's Pre-Installation Manual for exact information regarding fault expansions).

coordination of the magnetic shield design into existing or planned facilities, such as foundations/footings and other building

agents, and/or representatives as part of the MR Site Evaluation process. The final decision on the Design Objectives criteria and the

GE recommends that your final shield fabrication drawings be forwarded for review. Any deviations from the design specified in this drawing set must be brought to the attention of GE. Please allow up to two weeks for the review process.

Send drawings to your Project Manager as identified on the cover sheet.

ColumbiaUMC_20221448_run00_NewYorkNY_15IPM-vf.xlsm

SIGNA ARTIST | 1:100 | Rev A | Date | 30/JUL/2020 MRI-M221448-FIN-00.DWG

magnetic shielding design applies solely to the magnetic resonance suite configured as shown in this drawing set.

Contain the MRI-generated magnetic induction of 5 gauss within 11'-10" plan north of the magnet's isocenter.

- magnetic field and;
- Customer demonstrates to GE that the shielding material and fabrications are in strict accordance with GE's design and
- 4. Customer demonstrates to GE that the shielding is installed in strict accordance with GE's design, specifications and
- 5. Customer demonstrates to GE that (i) the System has not been and is not being subjected to improper or extraordinary use, (ii)

The Customer is responsible for the effect of the fringe fields produced by the System's magnet, and is responsible for the components.

The Design Objectives have been formulated by GE from information, both written and verbal, obtained from the Customer, his amount of magnetic field containment needed at this site rests solely with the Customer or his designated representatives. Please contact MR Siting/Shielding immediately with any corrections and/or additions related to this magnetic shield design.

A9 - Shield Requirements (1)

11/27

SHEET TITLE: GE MR-A9

SEAL:

DATE:06/18/2021 CON/REF No. CONTRACT No. **SCALE:** AS NOTED **PROJECT No.** 6109 CHECKED: CH DRAWN: KU

DATE

Architect of Record:

ARRAY-ARCHITECTS.COM

STRUCTURAL ENGINEER

LORING CONSULTING ENGINEERS, INC.

NEW YORK, NY. 10016

212-689-3110

CONSULTANTS:

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER

360 WEST 31ST STREET

NEW YORK, NY 10001 PHONE: 646.674.6100

OWNER:

PROJECT:

NEW MRI

TARRYTOWN

155 WHITE PLAINS ROAD

TARRYTOWN, NY 10591

KEY PLAN

COLUMBIA DOCTOR'S

W - WEST BUILDING E - EAST BUILDING

470 PARK AVE SOUTH, 11th FLOOR

SHEET NO. V-111.00

DESCRIPTION

REVISIONS/ISSUES

DWG OF 47

CD Submission 06/18/2021

This document is a copyright protected instrument of service, prosecuted under the law.

C COPYRIGHT - 2020 ARRAY-ARCHITECTS P.C

property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be

470 PARK AVE SOUTH, 11th FLOOR

NEW YORK, NY. 10016 212-689-3110

CONSULTANTS:

STRUCTURAL ENGINEER

ARRAY-ARCHITECTS.COM

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER

PHONE: 646.674.6100

LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001


OWNER:


COLUMBIA DOCTOR'S TARRYTOWN

PROJECT: NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591







SHEET TITLE: GE MRI-A10


SEAL:

CON/REF No. CONTRACT No. SCALE: AS NOTED PROJECT No. 6109 CHECKED: CH DRAWN: KU

DATE:06/16/21

SHEET NO. V-112.00

DWG OF 47 This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be prosecuted under the law.



CONSULTANTS:

STRUCTURAL ENGINEER

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER


LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001 PHONE: 646.674.6100


OWNER:

COLUMBIA DOCTOR'S TARRYTOWN

PROJECT: **NEW MRI**

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591





DATE **DESCRIPTION REVISIONS/ISSUES**

SHEET TITLE: GE MRI-A11

SEAL:

CON/REF No. CONTRACT No. **SCALE:** AS NOTED **PROJECT No.** 6109 CHECKED: CH DRAWN: KU

V-113.00

DWG OF 47

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be




ARMCO DI-MAX M36 (or material with equivalent magnetic properties) CR FP BARE (non coated) non-oriented electrical steel (2.25% maximum silicon content, 0.004% maximum carbon content, 0.65% maximum aluminum content).


- Material must meet all requirements of ASTM A677 for fully-processed non-oriented electrical steel.
- Material must meet or exceed magnetic properties as illustrated in figure 1 when tested per ASTM A596 as sheared with laminations in test specimen half parallel, half transverse to the rolling direction.
- Material must not be handled with electromagnetic equipment.
- Each sheet of material must lie flat under it's own weight. Material supplier must provide manufacturers certification with shipment.

- Magnetic properties of electrical steels are especially sensitive to stress, substantial reductions in magnetic properties can be caused by only a few
- hundred pounds per square inch. therefore, prevention of stresses due to elastic or plastic deformation during handling is critical.
- Refer to figure 2 for some examples of situations to avoid.
- Each 1/8" (3 mm) shield thickness consists of 5 layers of 24 gauge (0.0250 inch [0.64mm]) M36 steel.

(SILICON-011310)




Excessive stresses may be produced in electrical steel by improper arrangement of coils in storage, as indicated.



Unintentional introduction of excessive stress in steel during storage of coils.

Impproper handling of strip, sheets or long laminations can introduce stress that impair magnetic properties of the steel.



Coil set introduced in handling strip or long laminations.

FIGURE 2

Columbia University Medical Center SIGNA ARTIST

MRI-M221448-FIN-00.DWG

| 1:100 |Rev A|Date | 30/JUL/2020

A11 - Shield Requirements (3)

13/27

1. Magnetic shield ceilings, floors, walls, etc. Will consist of multiple layers to make of the design's required thickness. 2. Slotted or oversized bolt holes are recommended to assure adequate play for positioning plates tightly together. Refer to detail a.

must be assessed and documented by GE MR Siting.

non-magnetic stainless steel is acceptable.

Siting and Planning.

3. If multiple layers of steel are used, some lapping plates can be eliminated by staggering plate separations. Refer to detail B. This is the recommended method for multiple layers.

4. If layers are not staggered or a single layer of the required thickness is used, lapping plate must be used to cover plate separation. Lapping plate is required at each corner, plate separations perpendicular to the magnet bore's axis. This include all vertical seams in the sidewalls, all

cross seams in floor and ceiling, all horizontal and vertical seams in the front and rear surfaces and all corners where a front or rear surface meet any other surface. Refer to details C, D, E, F, and shield isometric. 5. Gaps in the shield which are perpendicular to the flux path are critical and should be kept to a minimum. Maximum separation between plates

SHIELD NOTES

MAGNETIC SHIELD CONSTRUCTION REQUIREMENTS

The actual performance of the magnetic shield depends upon strict adherence to the following requirements.

GENERAL

2. The magnetic shield design must be adhered to exactly. Any deviation from dimensional design, fabrication technique or material specification

3. Building and shield support structures should be outside of the steel plates to preclude detrimental effects on homogeneity. If it is absolutely

4. For all 1.5T and 3T MR systems, use of non-magnetic materials for concrete floor reinforcement is recommended, unless a homogeneity

necessary to have support structure on the inside of the shield, a homogeneity analysis must be performed and documented by GE MR siting.

analysis has been performed and documented by GE MR Siting. For older 1.5T systems with the Magnishield option, signa 0.5T systems, MR

The minimum finished floor to finished ceiling clearance in the magnet room is shown on sheet A2. Clearance dimensions do not include space

5. A secondary plate must be magnetically isolated from the primary shield. The use of non-magnetic structural members such as aluminum or

6. If the magnetic shield is installed inside an RF shield, both the shield and its support system must be isolated from ground by 1,000 ohms or

8. Locations of the penetration panel, window, and door openings are established on the final layout. Changes must be documented by GE MR

11. There must be a minimum of 8 inch [200mm] between the magnetic shield and the area where the magnetic field is being reduced. For example

if a hallway has a 5 gauss restriction goal, the magnetic shield should be positioned 8 inch [200mm] inside the hallway surface. Larger buffer

possible. If the parent wall structure is steel, then the magnetic shield needs to be decoupled from the steel studs by a minimum of 5/8 inch

contacts. Refer to the RF section of the pre-installation manual for all requirements. Magnetic decoupling media may also be needed if the

FABRICATION

12. If the magnetic shield is supported by the parent room structure, the parent room structure should be wood/concrete/aluminum whenever

13. If the magnetic shield is supported inside the RF structure, care must be taken to avoid loose metal to metal contacts or dissimilar metal

The design of the shield support system is the responsibility of the customer's architect, structural engineer, and shield vendor.

Max systems and Vectra systems, conventional reinforcing and structural members are acceptable.

required for the magnetic shield, RF shield, shield structure, and finished ceiling materials.

9. Penetration panel mounting details are specified in the appropriate GE site planning direction.

Consult GE MR Siting for guidance on any issues relating to magnetic shielding.

spaces are preferred but 8 inch [200mm] is the recommended minimum.

magnetic shield is attached directly to an RF shield that is primarily steel.

[16mm] gypboard, plywood, furring, or other non-magnetic media.

greater and meet all the GE RF room requirements as specified in the appropriate GE pre-installation manual.

6. Hole diameter of 8" [200mm] or less, when separated from another hole by a minimum of 3 times the diameter of the hole is allowed as needed. Holes for cryogen venting, HVAC, etc. Must be assessed and documented by GE MR Siting.

7. Angled joints do not need to have beveled mating edges. Refer to detail E.

8. When lapping plate material is required, it should extend at least 3" [76mm] to each side of the separation, and should have a total minimum thickness of 1/3 the thickness of the shield thickness. 9. Any separation between plates which exceed 1/8" [3mm] cannot be corrected by placing a lapping plate over the separation or using full seam

welds across the separation. Corrections for this situation must be assessed and documented by GE MR Siting. 10. Lapping plate can be attached by tack welding or bolting.

11. Shield plates are to be free of loose scale. Failure to remove scale can result in loose scale and steel particles becoming lodged in the magnet.

The preferred method for removal of scale is sandblasting. Removal of scale using wire brush is permitted. Shot blasting is prohibited. 12. Shield plates must be rigidly supported to prevent movement resulting from air pressure or other environmentally induced changes which can alter the magnet's homogeneity or systems performance.

13. Shield plate must be well anchored to the support structure by plug welding or bolting.

14. Full seam welding is not recommended for attaching shield plates to support structure or adjacent shield plates.

15. Minimum penetration full seam welds are allowed when the magnetic shield also serves as the RF shield.

Full seam welding is not allowed on silicon steel sheets. 17. Minimize the need for flame cutting of steel plate after annealing, as material's magnetic properties can be altered as a result of stress induced

18. Loose ferrous components become dangerous projectiles when the magnet is energized.

CD Submission 06/18/2021

DATE:06/16/21

SHEET NO.

STRUCTURAL NOTES

- All units that are wall mounted or wall supported are to be provided with supports where necessary. Wall supports are to be supplied and installed by the customer or his contractors.
- Dimensions are to finished surfaces of room.
- Certain MR procedures require an extremely stable environment to achieve high resolution image quality. Vibration is known to introduce field instabilities into the imaging system. The vibration effects on image quality can be minimized during the initial site planning of the mr suite by minimizing the vibration environment. See PROXIMITY LIMITS, PATIENT TABLE DOCK ANCHOR MOUNTING REQUIREMENTS AND VIBROACOUSTIC DAMPENING KIT details for additional information.
- Standard steel studs, nails, screws, conduit, piping, drains and other hardware are acceptable if properly secured. Any loose steel objects can be violently accelerated into the bore of the magnet. Careful thought should be given to the selection of light fixtures, cabinets, wall decorations, etc. To minimize this potential hazard. For safety, all removable items within the magnet room such as faucet handles, drain covers, switch box cover plates, light fixture components, mounting screws, etc. must be non-magnetic. If you have a specific question about material, bring it to the attention of your GE project manager of installations.
- Floor levelness refer to MAGNET ROOM FLOOR SPECIFICATIONS DETAIL, this floor levelness requirement is important for accurate patient table docking.
- Non-movable steel such as wall studs or hvac components will produce negligible effect on the active shield
- Customers contractor must provide all penetrations in post tension floors.

Columbia University Medical Center


- Customers contractor must provide and install any non-standard anchoring. Documents for standard anchoring methods are included with GE equipment drawings for geographic areas that require such documentation.
- Customers contractor must provide and install hardware for "through the floor" anchoring and/or any bracing under access floors. This contractor must also provide floor drilling that cannot be completed because of an obstruction encountered while drilling by the GE installer such as rebar etc.
- Customers contractor to provide and install appropriate supports for the storage of excess cables.
- It is the customer's responsibility to perform any floor or wall penetrations that may be required. The customer is also responsible for ensuring that no subsurface utilities (e.g., electrical or any other form of wiring, conduits, piping, duct work or structural supports (i.e. post tension cables or rebar)) will interfere or come in contact with subsurface penetration operations (e.g. drilling and installation of anchors/screws) performed during the installation process. To ensure worker safety, GE installers will perform surface penetration operations only after the customer's validation and completion of the "GE surface penetration

SIGNA ARTIST

VIBRATION SPECIFICATIONS

Excessive vibration can affect MR image quality. Vibration testing must be performed early in the site planning process to ensure vibration is minimized. Both steady state vibration (exhaust fans, air conditioners, pumps, etc.) and transient vibrations (traffic, pedestrians, door slamming, etc.) must be assessed. The Magnet cannot be directly isolated from vibration. Any vibration issue must be resolved at the source.

Transient vibration levels above the specified limits in the MR Site Vibration Test Guidelines must be analyzed. Any transient vibration that causes vibration to exceed the steady-state level must be mitigated.



VIBRATION TRANSMITTED THROUGH VIBROACOUSTIC MAT 0.225 0.025 1/3 Octave Frequency [Hz] MRI-M221448-FIN-00.DWG 14/27 Rev A Date 30/JUL/2020 C1 - Cover Sheet

OWNER:

COLUMBIA DOCTOR'S TARRYTOWN

Architect of Record:

ARRAY-ARCHITECTS.COM

NEW YORK, NY. 10016

STRUCTURAL ENGINEER

LORING CONSULTING ENGINEERS, INC.

212-689-3110

CONSULTANTS:

REUTHER + BOWEN

DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER


360 WEST 31ST STREET


NEW YORK, NY 10001 PHONE: 646.674.6100

470 PARK AVE SOUTH, 11th FLOOR

PROJECT: NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591





DATE **DESCRIPTION REVISIONS/ISSUES**

SHEET TITLE: GE MR-S1

SEAL:

DATE:06/16/21 CON/REF No. CONTRACT No. **SCALE:** AS NOTED PROJECT No. 6109 CHECKED: CH DRAWN: KU

SHEET NO. V-114.00

DWG OF 47 This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title

470 PARK AVE SOUTH, 11th FLOOR

NEW YORK, NY. 10016 212-689-3110

CONSULTANTS:

STRUCTURAL ENGINEER

ARRAY-ARCHITECTS.COM

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER


LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001 PHONE: 646.674.6100

OWNER:

COLUMBIA DOCTOR'S TARRYTOWN

PROJECT: NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591





SHEET TITLE:
GE MRI-S2

KEY PLAN

SEAL: **DATE:**06/16/21 CON/REF No. CONTRACT No. SCALE: AS NOTED PROJECT No. 6109 CHECKED: CH DRAWN: KU


SHEET NO.

V-115.00

C COPYRIGHT - 2020 ARRAY-ARCHITECTS P.C.

DWG OF 47

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be prosecuted under the law.



CONSULTANTS:

STRUCTURAL ENGINEER

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER


LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001 PHONE: 646.674.6100




COLUMBIA DOCTOR'S TARRYTOWN

PROJECT:
NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591







SHEET TITLE:
GE MRI-S3

SEAL:


CON/REF No.
CONTRACT No.
SCALE: AS NOTED
PROJECT No. 6109
CHECKED: CH
DRAWN: KU

DATE:06/18/2021

V-116.00

DWG OF 47

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be prosecuted under the law.



DESCRIPTION

Emergency exhaust vent - refer to magnet room vent requirements (position to be defined)

Cryogen vent (200mm [8"] O.D.)

ARRAY-ARCHITECTS.COM
470 PARK AVE SOUTH, 11th FLOOR
NEW YORK, NY. 10016
212-689-3110

CONSULTANTS:

STRUCTURAL ENGINEER

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER


PHONE: 646.674.6100


LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001

OWNER:
COLUMBIA DOCTOR'S
TARRYTOWN

PROJECT:
NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591





NO. DESCRIPTION DATE
REVISIONS/ISSUES

SHEET TITLE: GE MRI-M1

SEAL:

DATE:06/18/2021

CON/REF No.

CONTRACT No.

SCALE: AS NOTED

PROJECT No. 6109


CHECKED: CH DRAWN: KU

SHEET NO.

V-117.00

DWG OF 47

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be



MECHANICAL/PLUMBING NOTES

_

С

_Columbia_Doctors_Tarrytown__

:8:28 PM 12" = 1'-0" C:\Proje

X N

Α

TEMPERATURE AND HUMIDITY SPECIFICATIONS

IN-USE CONDITIONS

	MA	AGNET RO	OM	co	NTROL ROC	MC	EQUIPMENT ROOM				
		Range			Range		Range				
Temperature		15 to 21°C			15 to 32°C		15 to 32°C				
		59 to 69.8°F			59 to 89.6°F		59 to 89.6°F				
Temperature gradient		± 3°C/h			± 3°C/h		± 3°C/h				
		± 5°F/h			± 5°F/h		± 5°F/h				
Relative humidity (1)	1	30% to 60%	7.		30% to 70%		30% to 70%				
Humidity gradient	≤ 5%/h			,	≤ 5%/h			≤ 5%/h			
	Stand by	Average	Max	Stand by	Average	Max	Stand by	Average	Max		
System heat dissipation	1.01kW	1.8kW	3.15kW	7.	1.46kW		5.79kW	6.87kW	13.05kW		
	3450 btu	6142 btu	10748 btu		4947 btu		19769 btu	23225 btu	44523 btu		

NOTE

Maximum ambient temperature for the Equipment room at inlet is derated by 1°C per 300 m (984 ft) above 2000 m (6562 ft) (not to exceed 2600 m [8530 ft]).

AIR EXCHANGE

Columbia University Medical Center

According to local standards.

NOTE

In case of using air conditioning systems or chilled water piping that have a risk of water leakage it is recommended not to install it above electric equipment or to take measures to protect the equipment from dropping water.

HEAT DISSIPATION DETAILS

DESCRIPTION	BOOM	ID	LE	AVE	RAGE	MAX	
DESCRIPTION	ROOM	W	btu	W	btu	W	btu
Magnet (MAG) and Patient Table (PT)	Magnet	561	1915	1200	4095	2400	8189
Blower Box (MG6)	Magnet	450	1535	450	1535	450	1535
Penetration Panel Cabinet (PEN)	Magnet	0	0	150	512	300	1024
Penetration Panel Cabinet (PEN)	Equipment	1568	5349	1568	5349	3135	10697
Secondary Penetration Wall (SPW)	Magnet/Equipment		0.00	AV 19	0		
Main Disconnect Panel (MDP)	Equipment	132	450	132	450	264	901
Power, Gradient, RF Cabinet (PGR)	Equipment	2500	8530	3068	10470	6137	20940
Crycooler Compressor (CRY)	Equipment	500	1706	500	1706	500	1706
Heat Exchanger Cabinet (HEC)	Equipment	500	1706	500	1706	1000	3412
Magnet Monitor (MON)	Equipment	240	819	240	819	240	819
Operator Workspace equipment (OW)	Control	1450	4947	1450	4947	1450	4947
OPTIONS							
MR Elastography (MRE)	Equipment	141	480	141	480	141	480

SIGNA ARTIST

MAGNET ROOM VENTING REQUIREMENTS

HVAC VENT REQUIREMENTS

- HVAC vendor must comply with Magnet room temperature and humidity specifications and RF shielding specifications.
- RF Shield vendor must install open pipe or honeycomb HVAC waveguides.
- All serviceable parts in the Magnet room (e.g.: diffusers) must be non-magnetic.
- Waveguides must be nonmagnetic and electrically isolated.
- Incoming air must contain at least 5% air from outside the Magnet room (inside or outside the facility) to displace residual helium.

EMERGENCY VENT REQUIREMENT

- Exhaust vent system is supplied by the customer.
- All items within the RF enclosure must be non-magnetic.
- The exhaust vent system must be tested and operational before the magnet is installed.
- The exhaust intake vent must be located near the magnet cryogenic vent at the highest point on the finished or drop ceiling.
- The Magnet room exhaust fan and exhaust intake vent must have a capacity of at least 1200 CFM (34 m³/min) with a minimum of 12 room air exchanges per hour.
- The exhaust fan must be placed above RF shielding located outside 10 gauss (1mT) and with appropriate waveguide.
- The system must have a manual exhaust fan switch near the Operator Workspace and in the Magnet room near the door (the switches must be connected in parallel).
- All system components must be accessible for customer inspection, cleaning and maintenance


PRESSURE VENT REQUIREMENT

MRI-M221448-FIN-00.DWG

- A pressure equalizing vent is required in the magnet room ceiling or in the wall, at the highest point possible.
- The vent minimum size must be (610 mm x 610 mm [24 in x 24 in]) or equivalent.
- The pressure equalization vent must be located so any Helium gas is not vented into occupied areas.

Rev AlDate 30/JUL/2020

Note: Location may affect acoustic noise transmission into occupied spaces.



ARRAY-ARCHITECTS.COM
470 PARK AVE SOUTH, 11th FLOOR
NEW YORK, NY. 10016

Architect of Record:

CONSULTANTS:

212-689-3110

STRUCTURAL ENGINEER

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER


LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001 PHONE: 646.674.6100

OWNER:

COLUMBIA DOCTOR'S TARRYTOWN

PROJECT: NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591



NO. DESCRIPTION DATE
REVISIONS/ISSUES

SHEET TITLE:
GE MR-M2

SEAL:

DATE:06/18/2021
CON/REF No.
CONTRACT No.
SCALE: AS NOTED
PROJECT No. 6109
CHECKED: CH
DRAWN: KU

V-118.00

DWG OF 47

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be

C COPYRIGHT - 2020 ARRAY-ARCHITECTS P.C.

18/27

C1 - Cover Sheet



CONSULTANTS:

STRUCTURAL ENGINEER

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER


LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001 PHONE: 646.674.6100


OWNER:

COLUMBIA DOCTOR'S TARRYTOWN

PROJECT: NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591







SHEET TITLE: GE MRI-M3

SEAL:


DATE:06/18/2021 CON/REF No. CONTRACT No. SCALE: AS NOTED PROJECT No. 6109 CHECKED: CH DRAWN: KU

SHEET NO.

V-119.00

DWG OF 47

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be prosecuted under the law.



CONSULTANTS:

STRUCTURAL ENGINEER

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER


LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001 PHONE: 646.674.6100

OWNER:

COLUMBIA DOCTOR'S TARRYTOWN

PROJECT:
NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591



NO. DESCRIPTION DATE
REVISIONS/ISSUES

SHEET TITLE:
GE MRI-M4

KEY PLAN

SEAL:

CON/REF No.
CONTRACT No.
SCALE: AS NOTED
PROJECT No. 6109
CHECKED: CH
DRAWN: KU

DATE:06/18/2021

V-120.00


DWG OF 47

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without

CD Submission 06/18/2021

20/27


C1 - Cover Sheet



SIGNA ARTIST

MRI-M221448-FIN-00.DWG

Columbia University Medical Center



Rev A Date 30/JUL/2020

470 PARK AVE SOUTH, 11th FLOOR

ARRAY-ARCHITECTS.COM

NEW YORK, NY. 10016

212-689-3110

CONSULTANTS:

STRUCTURAL ENGINEER

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER


LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001 PHONE: 646.674.6100


OWNER:

COLUMBIA DOCTOR'S TARRYTOWN

PROJECT: NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591





NO.	DESCRIPTION	DATE

SHEET TITLE: GE MRI-M5

SEAL:

CON/REF No. CONTRACT No. SCALE: AS NOTED PROJECT No. 6109 CHECKED: CH DRAWN: KU

SHEET NO. V-121.00

DWG OF 47

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be prosecuted under the law.

C COPYRIGHT - 2020 ARRAY-ARCHITECTS P.C.

MAGNET CRYOGENIC VENT SYSTEM PRESSURE DROP MATRIX

Outer dia. of pipe (D)	of sys	tem conen com gnet	nt Pressure drop for straight pipe					Std sweep 90° elbow				90° miter		
	ft	m	psi/ft	kPa/m	psi	kPa	psi	kPa	psi	kPa	psi	kPa	psi	kPa
	10	3.048	0.146	3.311	0.804	5.544	0.536	3.696	1.507	10.394	1.005	6.93	3.014	20.789
6 in.	20	6.096	0.253	5.715	1.356	9.355	0.904	6.237	2.543	17.54	1.696	11.694	5.087	35.081
(Inside	30	9.144	0.374	8.451	1.845	12.727	1.23	8.485	3.46	23.864	2.307	15.909	6.92	47.727
dia. ≥144mm)	40	12.192	0.473	10.699	2.278	15.708	1.518	10.472	4.271	29.453	2.847	19.635	8.541	58.906
	50	15.24	0.554	12.534	2.66	18.342	1.773	12.228	4.987	34.392	3.325	22.928	9.974	68.783
	60	18.288	0.62	14.019	2.997	20.668	1.998	13.779	5.619	38.753	3.746	25.835	11.238	77.506
	10	3.048	0.025	0.564	0.188	1.294	0.125	0.862	0.352	2.426	0.234	1.617	0.703	4.851
	20	6.096	0.043	0.97	0.313	2.158	0.209	1.439	0.587	4.046	0.391	2.697	1.173	8.092
	30	9.144	0.064	1.45	0.427	2.944	0.285	1.963	0.8	5.52	0.534	3.68	1.601	11.04
	40	12.192	0.082	1.862	0.53	3.658	0.354	2.439	0.995	6.859	0.663	4.573	1.989	13.718
8 in.	50	15.24	0.098	2.215	0.624	4.307	0.416	2.871	1.171	8.075	0.781	5.383	2.342	16.15
(Inside dia.	60	18.288	0.111	2.516	0.71	4.895	0.473	3.263	1.331	9.179	0.887	6.119	2.662	18.357
≥195mm)	80	24.384	0.132	2.987	0.857	5.914	0.572	3.942	1.608	11.088	1.072	7.392	3.216	22.176
	100	30.48	0.147	3.318	0.979	6.752	0.653	4.501	1.836	12.659	1.224	8.439	3.671	25.318
	120	36.576	0.157	3.545	1.079	7.44	0.719	4.96	2.023	13.95	1.3490	9.3	4.046	27.901
	140	42.672	0.163	3.693	1.161	8.006	0.774	5.338	2.177	15.012	1.451	10.008	4.353	30.024
	160	48.768	0.167	3.784	1.228	8.471	0.819	5.648	2.303	15.884	1.535	10.589	4.606	31.768
	180	54.864	0.169	3.833	0.099	8.853	0.856	5.902	2.407	16.6	1.605	11.067	4.814	33.2
	40	6.096 12.192	0.011	0.241	0.168	0.683 1.16	0.066	0.455	0.186	1.28 2.175	0.124	0.854 1.45	0.371	2.561 4.351
	60	18.288	0.021	0.468	0.227	1.568	0.112	1.045	0.426	2.173	0.21	1.96	0.853	5.88
	80	24.384	0.025	0.781	0.278	1.916	0.132	1.277	0.521	3.592	0.347	2.395	1.042	7.184
10 in.	100	30.48	0.039	0.884	0.321	2.212	0.183	1.474	0.601	4.147	0.401	2.765	1.203	8.294
(Inside dia.	120	36.576	0.042	0.96	0.357	2.464	0.238	1.642	0.67	4.619	0.447	3.08	1.34	9.239
≥247mm)	140	42.672	0.042	1.016	0.388	2.678	0.259	1.785	0.728	5.021	0.485	3.348	1.456	10.043
	160	48.768	0.047	1.056	0.415	2.86	0.276	1.907	0.778	5.363	0.518	3.575	1.555	10.726
	180	54.864	0.048	1.083	0.437	3.015	0.291	2.01	0.82	5.653	0.547	3.769	1.64	11.307
	200	60.99	0.049	1.101	0.456	3.147	0.304	2.098	0.856	5.9	0.57	3.934	1.711	11.801
	20	6.096	0.004	0.08	0.04	0.277	0.027	0.184	0.075	0.519	0.05	0.346	0.15	1.037
	40	12.192	0.007	0.157	0.068	0.47	0.045	0.313	0.128	0.88	0.085	0.587	0.255	1.761
	60	18.288	0.01	0.22	0.093	0.638	0.062	0.425	0.174	1.197	0.116	0.798	0.347	2.393
24.9 (2)	80	24.384	0.012	0.269	0.114	0.786	0.076	0.524	0.214	1.473	0.142	0.982	0.427	2.946
12 in. (Inside	100	30.48	0.014	0.309	0.133	0.914	0.088	0.609	0.248	1.714	0.166	1.142	0.497	3.427
dia.	120	36.576	0.015	0.34	0.149	1.026	0.099	0.684	0.279	1.923	0.186	1.282	0.558	3.847
≥298mm)	140	42.672	0.016	0.364	0.163	1.123	0.109	0.749	0.305	2.106	0.204	1.404	0.611	4.212
	160	48.768	0.017	0.382	0.175	1.208	0.117	0.805	0.328	2.265	0.219	1.51	0.657	4.529
	180	54.864	0.018	0.396	0.186	1.282	0.124	0.854	0.348	2.403	0.232	1.602	0.697	4.806
	200	60.96	0.018	0.406	0.195	1.346	0.13	0.897	0.366	2.523	0.244	1.682	0.732	5.046

Outer dia. of pipe (D)	of vent system componen t from magnet		drop for straight pipe		Std sweep 45° elbow 45° el		200 M					90°	miter	
	ft	m	psi/ft	kPa/m	psi	kPa	psi	kPa	psi	kPa	psi	kPa	psi	kPa
	20	6.096	0.001	0.032	0.019	0.13	0.013	0.086	0.035	0.243	0.024	0.162	0.071	0.486
	40	12.192	0.003	0.063	0.032	0.219	0.021	0.146	0.06	0.411	0.04	0.274	0.119	0.823
	60	18.288	0.004	0.088	0.043	0.299	0.029	0.2	0.081	0.561	0.054	0.374	0.163	1.122
14 in. (Inside dia. ≥348mm)	80	24.384	0.005	0.11	0.054	0.37	0.036	0.247	0.101	0.694	0.067	0.463	0.201	1.389
	100	30.48	0.006	0.127	0.063	0.433	0.042	0.289	0.118	0.812	0.079	0.542	0.236	1.625
	120	36.576	0.006	0.141	0.071	0.489	0.047	0.326	0.133	0.917	0.089	0.611	0.266	1.834
	140	42.672	0.007	0.152	0.078	0.539	0.052	0.359	0.146	1.01	0.098	0.673	0.293	2.02
	160	48.768	0.007	0.161	0.084	0.582	0.056	0.388	0.158	1.092	0.106	0.728	0.317	2.184
ì	180	54.864	0.007	0.169	0.090	0.621	0.06	0.414	0.169	1.164	0.113	0.776	0.338	2.329
8	200	60.96	0.008	0.174	0.095	0.655	0.063	0.437	0.178	1.229	0.119	0.819	0.356	2.457
38	20	6.096	0.001	0.014	0.01	0.068	0.007	0.045	0.018	0.127	0.012	0.084	0.037	0.253
5	40	12.192	0.001	0.028	0.017	0.114	0.011	0.076	0.0310	0.2130	0.021	0.142	0.062	0.427
3	60	18.288	0.002	0.04	0.023	0.156	0.015	0.104	0.042	0.292	0.028	0.195	0.085	0.584
16 in.	80	24.384	0.002	0.05	0.028	0.193	0.019	0.129	0.053	0.362	0.035	0.242	0.105	0.725
(Inside	100	30.48	0.003	0.059	0.033	0.227	0.022	0.151	0.062	0.426	0.041	0.284	0.124	0.852
dia. ≥399mm)	120	36.576	0.003	0.066	0.037	0.258	0.025	0.172	0.07	0.483	0.047	0.322	0.14	0.966
233311111)	140	42.6720	0.003	0.072	0.041	0.285	0.028	0.19	0.077	0.534	0.052	0.356	0.155	1.069
8	160	48.768	0.003	0.076	0.045	0.309	0.03	0.206	0.084	0.58	0.056	0.387	0.168	1.16
	180	54.864	0.004	0.08	0.048	0.331	0.032	0.221	0.09	0.621	0.06	0.414	0.18	1.243
	200	60.96	0.004	0.083	0.051	0.351	0.034	0.234	0.095	0.658	0.064	0.439	0.191	1.316

- Elbows with angles greater than 90 deg must not be used
- Data in Table is based on the following facts and assumptions:
- Initial flow conditions at magnet interface
- EM energy (13MJ) is dumped to He during quench and rises He temperature to 10 Kelvin Gas temperature starting at 10 Kelvin and increase with length determined by thermal energy balance
- 90% He is assumed to be evacuated within 30 sec. None left after quench.
- Absolute roughness is assumed to be 0.25 mm.
- R/D = 1.0 for standard sweep elbows, R/D = 1.5 for long sweep elbows where D = outer diameter of pipe;
- R = radius of bend

The total pressure drop of the entire cryogenic vent system must be less than 20 psi (137.9 kPa). The calculation starts at the magnet vent interface and ends at the termination point outside the building.

Columbia University Medical Center

SIGNA ARTIST

MRI-M221448-FIN-00.DWG

|Rev A|Date 30/JUL/2020 |

C1 - Cover Sheet

| 21/27

DATE:06/18/2021

 All lighting fixtures and associated components must meet all RF shielded room and RF grounding requirements (e.g., track lighting is not recommended due to possible RF noise).

- All removable lighting fixtures and associated components must be non-magnetic.
- All lighting must use direct current (the DC must have less than 5% ripple).
- 300 lux must be provided at the front of the magnet for patient access and above the magnet for servicing.

LIGHTING REQUIREMENTS

- Fluorescent lighting must not be used in the magnet room.
- Lighting must be adjusted using a discrete switch or a variable DC lighting controller.
- SCR dimmers or rheostats must not be used.
- DC LED lighting may be used if the DC power converter and RF sources are all located outside the magnet room RF shield.

NOTE: LED lighting could cause image quality issues due to RF interference. Make sure a MR-compatible LED lighting solution is chosen.

- Battery chargers (e.g., used for emergency lighting) must be located outside the magnet room.
- Short filament length bulbs are recommended.
- Linear lamps are not recommended due to the high burnout rate.

CONNECTIVITY REQUIREMENTS

Broadband Connections are necessary during the installation process and going forward to ensure full support from the Engineering Teams for the customers system. Maximum performance and availability for the customers system is maintained and closely monitored during the lifetime of the system. Proactive and reactive maintenance is available utilising the wide range of digital tools using the connectivity solutions listed below:

- Site-to-Site VPN/GE Solution
- Site-to-Site VPN/Customer Solution
- Connection through Dedicated Service Network
- Internet Access connectivity for InSite 2.0

The requirements for these connectivity solutions are explained in the broadband solutions catalogue (separate document).

ELECTRICAL NOTES

- All wires specified shall be copper stranded, flexible, thermo-plastic, color coded, cut 10 foot long at outlet boxes, duct termination points or stubbed conduit ends. All conductors, power, signal and ground, must be run in a conduit or duct system. Electrical contractor shall ring out and tag all wires at both ends. Wire runs must be continuous copper stranded and free from splices.
- 1.1. Aluminum or solid wires are not allowed.
- Wire sizes given are for use of equipment. Larger sizes may be required by local codes.
- It is recommended that all wires be color coded, as required in accordance with national and local electrical codes.
- Conduit sizes shall be verified by the architect, electrical engineer or contractor, in accordance with local or national codes.
- Convenience outlets are not illustrated. Their number and location are to be specified by others. Locate at least one convenience outlet close to the system control, the power distribution unit and one on each wall of the procedure room. Use hospital approved outlet or equivalent.
- 6. General room illumination is not illustrated. Caution should be taken to avoid excessive heat from overhead spotlights. Damage can occur to ceiling mounting components and wiring if high wattage bulbs are used. Recommend low wattage bulbs no higher than 75 watts and use dimmer controls (except MR). Do not mount lights directly above areas where ceiling mounted accessories will be parked.
- Routing of cable ductwork, conduits, etc., must run direct as possible otherwise may result in the need for greater than standard cable lengths (refer to the interconnection diagram for maximum usable lengths point to point)
- 8. Conduit turns to have large, sweeping bends with minimum radius in accordance with national and local electrical codes.
- 9. A special grounding system is required in all procedure rooms by some national and local codes. It is recommended in areas where patients might be examined or treated under present, future, or emergency conditions. Consult the governing electrical code and confer with appropriate customer administrative personnel to determine the areas requiring this type of grounding system.
- 10. The maximum point to point distances illustrated on this drawing must not be exceeded.
- 11. Physical connection of primary power to GE equipment is to be made by customers electrical contractor with the supervision of a GE representative. The GE representative would be required to identify the physical connection location, and insure proper handling of GE equipment.
- 12. GEHC conducts power audits to verify quality of power being delivered to the system. The customer's electrical contractor is required to be available to support this activity.
- All junction boxes, conduit, duct, duct dividers, switches, circuit breakers, cable tray, etc., are to be supplied and installed by customers electrical contractor.
- Conduit and duct runs shall have sweep radius bends
- Conduits and duct above ceiling or below finished floor must be installed as near to ceiling or floor as possible to reduce run length.
- Ceiling mounted junction boxes illustrated on this plan must be installed flush with finished ceiling.
- All ductwork must meet the following requirements:
- 1.Ductwork shall be metal with dividers and have removable, accessible covers.
- 2.Ductwork shall be certified/rated for electrical power purposes.
- 3. Ductwork shall be electrically and mechanically bonded together in an approved manner.
- 4.PVC as a substitute must be used in accordance with all local and national codes. All openings in raceway and access flooring are to be cut out and finished off with grommet material by the
- General contractor to insert pull cords for all cable run conduits between the equipment room and the operators control room.
- 10 foot pigtails at all junction points.

customers contractor.

 Grounding is critical to equipment function and patient safety. Site must conform to wiring specifications shown on this plan.

| 22/27 Columbia University Medical Center MRI-M221448-FIN-00.DWG Rev A Date 30/JUL/2020 SIGNA ARTIST C1 - Cover Sheet

REUTHER + BOWEN

ARRAY-ARCHITECTS.COM

STRUCTURAL ENGINEER

NEW YORK, NY. 10016

212-689-3110

CONSULTANTS:

470 PARK AVE SOUTH, 11th FLOOR

Architect of Record:

MEP ENGINEER

DUNMORE, PA 18512 PHONE: 570.496.7020


LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001 PHONE: 646.674.6100

OWNER:

COLUMBIA DOCTOR'S TARRYTOWN

PROJECT: NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591



KEY PLAN

DESCRIPTION DATE **REVISIONS/ISSUES**

SHEET TITLE: GE MR-E1

SEAL:

DATE:06/18/2021 CON/REF No. CONTRACT No. **SCALE:** AS NOTED **PROJECT No.** 6109 CHECKED: CH DRAWN: KU

SHEET NO. V-122.00

DWG OF 47

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be

C COPYRIGHT - 2020 ARRAY-ARCHITECTS P.C

CD Submission 06/18/2021

Outlet Legend for GE Equipment

System emergency off (SEO), (recommended height 1.2m [48"] above floor)

Ooor interlock switch

DESCRIPTION

(CONTRACTOR SUPPLIED & INSTALLED)

ITEM

Cable ladder 450mm x 150mm [18" x 6"]

Cable ladder 450mm x 150mm [18" x 6"] for gradient cables

Additional Conduit Runs

(Contractor Supplied and Installed)

Power, Gradient, RF cabinet

From

Qty Size (in) Size (mm)

as Req'd

ARRAY-ARCHITECTS.COM 470 PARK AVE SOUTH, 11th FLOOR **NEW YORK, NY. 10016** 212-689-3110

CONSULTANTS:

STRUCTURAL ENGINEER

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER


LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001 PHONE: 646.674.6100

OWNER:

COLUMBIA DOCTOR'S TARRYTOWN

PROJECT: NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591



W - WEST BUILDING E - EAST BUILDING **KEY PLAN**

DESCRIPTION REVISIONS/ISSUES


SHEET TITLE:
GE MR-E2


SEAL:

DATE:06/18/2021 CON/REF No. CONTRACT No. SCALE: AS NOTED PROJECT No. 6109 CHECKED: CH DRAWN: KU

SHEET NO. V-123.00

DWG OF 47





470 PARK AVE SOUTH, 11th FLOOR

CONSULTANTS:

212-689-3110

STRUCTURAL ENGINEER

ARRAY-ARCHITECTS.COM

NEW YORK, NY. 10016

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER


LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001 PHONE: 646.674.6100


OWNER:

COLUMBIA DOCTOR'S TARRYTOWN

PROJECT: NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591







SHEET TITLE:
GE MRI-E3


SEAL: **DATE:**06/18/2021 CON/REF No. CONTRACT No. SCALE: AS NOTED PROJECT No. 6109 CHECKED: CH

DRAWN: KU

SHEET NO. V-124.00

DWG OF 47

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be prosecuted under the law.



CONSULTANTS:

STRUCTURAL ENGINEER

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER

LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001 PHONE: 646.674.6100


OWNER:

COLUMBIA DOCTOR'S TARRYTOWN

PROJECT: NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591







SHEET TITLE: GE MR-E4


SEAL:

DATE:06/18/2021 CON/REF No. CONTRACT No. SCALE: AS NOTED PROJECT No. 6109 CHECKED: CH DRAWN: KU

SHEET NO. V-125.00

DWG OF 47

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be prosecuted under the law.



POWER REQUIREMENTS SPECIFICATIONS OF MAIN POWER INPUT 380/400/415/480V ±10%, THREE-PHASE + N + G POWER SUPPLY **FREQUENCIES** 50/60Hz ± 3Hz

123kVA

< 17kVA

99kVA

- Power input must be separated from any others which may generate transients (elevators, air conditioning, radiology rooms equipped with
- high speed film changers...).

POWER FACTOR

INSTALLED LOAD

STAND-BY POWER

 Total harmonic distortion less than 2.5%. Phase imbalance must not exceed 2%

MAXIMUM INPUT POWER (5 sec MAX)

 Lock-out/Tag-out: The Main Disconnect Panel (MDP) shall provide an external single point lock-out/tag-out feature for the entire system and a means to externally lock-out/tag-out each output breaker independently. Each lock-out/tag-out feature shall accommodate a standard sized lock hasp.

SPECIFICATIONS OF BACK-UP POWER SUPPLY

MAGNET MONITOR REQUIRES A 110/220 VAC, 50/60 HZ, 2.0 A FACILITY SUPPLIED OUTLET. POWER AT THE **OUTLET MUST BE CONTINUOUSLY AVAILABLE.**

FOR CRYOCOOLER COMPRESSOR	
POWER INPUT	380/400/415/480V, THREE-PHASE + G
POWER REQUIREMENT	MIN 9kVA
POWER CONSUMPTION	MAX 7.2kW / STEADY STATE 6.5kW at 50Hz MAX 8.3kW / STEADY STATE 7.5kW at 60Hz
FREQUENCY	50/60Hz ± 3Hz

CABLES

- Power and cable installation must comply with the distribution diagram.
- Size of the Main power input cable is determined by the customer, taking its length and admissible voltage drops into consideration.
- All cables must be isolated and flexible, cable color codes must comply with standards for electrical installation.
- The cables from signaling and remote control (Y, Emergency Off Buttons, L...) will go to Main Panel with a pigtail length of 1.5m [60in], and will be connected during installation.
- Each conductor will be identified and isolated (screw connector).

GROUND SYSTEM


The equipotential link will be by means of an equipotential bar.

Columbia University Medical Center

- The grounding point of MDP is directly connected to the building's ground by an isolated copper cable.
- The impedance of the earth bar should be less than or equal to 2 ohms.

Direct feed from facility to MR system	MR system fed by dedicated fa	acility distribution transformer
Pro	erequisite Conditions	
480V 3-phase		
7,900 kVA (at source of feeder to MDP)	8,325 kVA (at input to distribution transformer)	
460V	475V (transformer secondary tapped accordingly)	
Feeder and 1	Transformer Recommendations	
N/A	Size: 225 kVA Impedance (Z): ≤5% K-Factor: ≥ K=20 200A overcurrent protection on secondary*	Size: 225 kVA Impedance (Z): ≤4% K-Factor: ≥ K=20 200A overcurrent protection on secondary*
280 ft	150 ft	240 ft
3/0 AWG Cu	3/0 AWG Cu	
6 AWG Cu (equipment grounding conductor)	4 AWG Cu (supply side bonding jumper) 6 AWG Cu (equipment grounding conductor)	
6 AWG Cu (bonding conductor)	6 AWG Cu (bonding conductor)	
	7,900 kVA (at source of feeder to MDP) 460V Feeder and 1 N/A 280 ft 3/0 AWG Cu 6 AWG Cu (equipment grounding conductor) 6 AWG Cu (bonding	Prerequisite Conditions 480V 3-phase 7,900 kVA (at source of feeder to MDP) 460V 475V (transformer secondary* N/A Size: 225 kVA Impedance (Z): ≤5% K-Factor: ≥ K=20 200A overcurrent protection on secondary* 280 ft 3/0 AWG Cu 6 AWG Cu (equipment grounding conductor) 6 AWG Cu (bonding 480V 3-phase 8,32: (at input to distrib 475V (transformer secondary* 150 ft 3/5 AWG Cu (supply secondary) 4 AWG Cu (supply secondary) 6 AWG Cu (equipment grounding conductor) 6 AWG Cu (bonding)

SIGNA ARTIST



EMERGENCY OFF LOCKS OUT ALL CONTRACTORS.

GE MDP IS UL AND cUL LABELED.

MRI-M221448-FIN-00.DWG

ALL CIRCUITS REQUIRE GROUND WIRES.

GE MDP SHORT CIRCUIT CURRENT RATING IS 25,000 AMPERES AT 480 VAC.

THE WIRE SIZE FOR THE EMERGENCY-OFF CIRCUIT IS 12-22 AWG CUSTOMER SUPPLIED

Rev AlDate 30/JUL/2020

Architect of Record:

ARRAY-ARCHITECTS.COM 470 PARK AVE SOUTH, 11th FLOOR **NEW YORK, NY. 10016** 212-689-3110

CONSULTANTS:

STRUCTURAL ENGINEER

REUTHER + BOWEN 326 WARD STREET DUNMORE, PA 18512 PHONE: 570.496.7020

MEP ENGINEER


LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001 PHONE: 646.674.6100

OWNER:

COLUMBIA DOCTOR'S TARRYTOWN

PROJECT: NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591



DATE **DESCRIPTION REVISIONS/ISSUES**

SHEET TITLE: GE MRI-E5

SEAL:

| 26/27

C1 - Cover Sheet

DATE:06/18/2021 CON/REF No. CONTRACT No. **SCALE:** AS NOTED PROJECT No. 6109 CHECKED: CH DRAWN: KU

SHEET NO. V-126.00

DWG OF 47

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be

CD Submission 06/18/2021

INTERCONNECTIONS

[94.5']

29.00m

S: 9.8m [32.2'] L: 16.8m [55.1'] 470 PARK AVE SOUTH, 11th FLOOR

212-689-3110

ARRAY-ARCHITECTS.COM

NEW YORK, NY. 10016


LORING CONSULTING ENGINEERS, INC. 360 WEST 31ST STREET NEW YORK, NY 10001 PHONE: 646.674.6100


OWNER:

COLUMBIA DOCTOR'S TARRYTOWN

PROJECT: NEW MRI

155 WHITE PLAINS ROAD TARRYTOWN, NY 10591







SHEET TITLE: GE MR-E6

SEAL:

CON/REF No. CONTRACT No. SCALE: AS NOTED PROJECT No. 6109 CHECKED: CH DRAWN: KU

DATE:06/18/2021

SHEET NO. V-127.00

DWG OF 47

This document is a copyright protected instrument of service, property of Array Architects and licensed for use in the title project only. Reproduction or use of this document without written permission of Array Architects is illegal and will be

